


Speciality Products



Unique procedure kits have been developed for Trans Oesophageal Echocardiography and Laparoscopic Ultrasound procedures.

For echo-endoscopes a complete line of water stand-off balloons is available.



ULTRACOVER®

Professional covers for ultrasound transducers

As you can see, ULTRACOVER® offers a solution for every application. If you need a special cover, do not hesitate to contact us. Together we will find the proper solution, one that meets the requirements of your specific application.

Microtek Medical is your trusted source of innovative product solutions for risk reduction and protection.

ULTRACOVER® offers the most complete range of disposable covers, products which facilitate the safe use of ultrasound probes.



P.O. Box 234 7200 AE Zutphen The Netherlands
Gerritsenweg 7 7202 BP Zutphen The Netherlands
Tel. +31 (0)575 599 200 Fax +31 (0)575 599 299

www.microtekmed.com

www.microtekmed.com



Covers for all branches of medicine

Whatever the branch of medicine, wherever ultrasound is used, ULTRACOVER® offers the proper solution. ULTRACOVER® is used in:

- Cardiology
- Gynecology
- Urology
- Radiology
- Intraoperative applications
- Gastroenterology



Our Endocavity line of covers offers a choice of materials to meet your performance criteria. Our latex products are dipped and rolled using medical-grade material. Our latex free IsoSilkí covers, for those with latex sensitivity, are folded to ease gel filling and application to the probe.

These covers are used in Transvaginal Imaging, Transrectal and Prostate Imaging and Biopsy.

This line of General Purpose Probe Covers, featuring IsoSilkí, is designed to enhance your general imaging needs when disease and contamination control is a high priority. All covers are sterile, telescopically folded, contain sterile gel and medical-grade tape strips, and latex free rubber bands to customize the fit.

These covers are used in Biopsy, Trauma-emergency room procedures, Patients in isolation, Catheter and drainage sites and to maintain sterility and cross-contamination.



Our Intraoperative Probe Covers incorporate a patented process allowing us to attach soft and supple, special-shaped probe tips to our latex free IsoSilkí fabric for full cable coverage. Special tips have been designed to provide better performance and manipulation. The tips are dipped for a perfect shape in our latex-free dipping material, or in natural latex. All surgical probe drapes are sterile, telescopically folded, contain sterile gel and medical grade tape strips, and latex-free bands to customize fit.

These covers are used in Abdominal Surgery, Cardiac Surgery and Vascular Surgery.



Our Surgical/Intraoperative Covers have been designed specifically for draping phased arrays, and curved linear arrays during surgical procedures. These drapes are long enough to cover the probe head and cable, all the way to the system connector, to guard against compromising the sterile field. The drapes are constructed using our latex-free IsoSilkí material. All surgical/intraoperative covers include latex-free bands, tape strips and sterile gel.

The Surgical/Intraoperative covers are used in Biopsy, Abdominal Surgery, Special Procedures, Cardiac Surgery and Vascular Surgery.