

Drains

Silicon Drains



Part.no.	Description	Qty/cs
Silicon Drains for general surgical use f.i. abdominal surgery		
Length 50 cm, 5 eyes with x-ray striping		
405108	Silicon Drains Ch. 8	25
405110	Silicon Drains Ch. 10	25
405112	Silicon Drains Ch. 12	25
405115	Silicon Drains Ch. 15	25
405118	Silicon Drains Ch. 18	25
405121	Silicon Drains Ch. 21	25
405124	Silicon Drains Ch. 24	25
405127	Silicon Drains Ch. 27	25
405130	Silicon Drains Ch. 30	25
405136	Silicon Drains Ch. 36	25

Silicon T-Drains



Part.no.	Description	Qty/cs
For f.i. galway drainage		
604003	Silicon T-Drains 3 mm	25
604004	Silicon T-Drains 4 mm	25
604005	Silicon T-Drains 5 mm	25
604006	Silicon T-Drains 6 mm	25
604007	Silicon T-Drains 7 mm	25
604008	Silicon T-Drains 8 mm	25

Flat Drains



Part.no.	Description	Qty/cs
Flat Drains type "Jackson Pratt" of 100% medical silicon, x-ray positive, double-layer peel-pack, sterile		
851204	Flat Drain, Width 4, Length 65 cm, 1-row eyes	50
851207	Flat Drain, Width 7, Length 125 cm, 2-rows eyes	50
851210	Flat Drain, Width 10, Length 105 cm, 2-rows eyes	50
Suction balloon		
851120	Silicon Bulb	25

Ecolab Healthcare
Clean. Safe. Healthy.

Microtek Medical B.V.
A Division of Ecolab
Hekkehorst 24
NL-7207 BN Zutphen
Tel.: +31 (0) 575 599200

Ecolab Europe GmbH
Richtistr. 7
CH-8304 Wallisellen
www.ecolab.com



Drains

Silicon Capillary-Drains



Part.no.	Description	Qty/cs
For general surgery f.i. hand surgery and gynaecology		
Thanks to the grooved internal wall, silicon capillary drains have an excellent capillary effect		
421406	Silicon Capillary-Drain, length 30 cm, width 6 mm with x-ray	25
421408	Silicon Capillary-Drain, length 30 cm, width 8 mm with x-ray	25
421410	Silicon Capillary-Drain, length 30 cm, width 10 mm with x-ray	25
421412	Silicon Capillary-Drain, length 30 cm, width 12 mm with x-ray	25
421612	Silicon Capillary-Drain, length 50 cm, width 12 mm with x-ray	25

Silicon Thorax-Drains



Part.no.	Description	Qty/cs
Specially designed for thoracic and heart surgery.		
These Drains can be connected perfectly with our Thoracic Drainage Systems		
401118	Silicon Thorax-Drains Ch.18, 50 cm with 8 eyes and 12 Perf. Cm and x-ray contrast	20
401124	Silicon Thorax-Drains Ch.24, 50 cm with 10 eyes and 17 Perf. Cm and x-ray contrast	20
401127	Silicon Thorax-Drains Ch.27, 60 cm with 10 eyes and 20 Perf. Cm and x-ray contrast	20
401130	Silicon Thorax-Drains Ch.30, 65 cm with 12 eyes and 24 Perf. Cm and x-ray contrast	20
401136	Silicon Thorax-Drains Ch.36, 70 cm with 14 eyes and 27 Perf. Cm and x-ray contrast	20

PVC Thorax-drains



Part.no.	Description	Qty/cs
Length 50 cm, 40 cm perforation (6 eyes) with x-ray contrast and atraumatic tip		
Deep marker every 2 cm		
401084	PVC Thorax-drain Ch.24, straight, with adapter	10
401094	PVC Thorax-drain Ch. 24, curved, trocar compliant	10
401088	PVC Thorax-drain Ch.28, straight, with adapter	10
401068	PVC Thorax-drain Ch. 28, straight, trocar compliant	10
401058	PVC Thorax-drain Ch.28, curved, with adapter	10
401098	PVC Thorax-drain Ch. 28, curved, trocar compliant	10
401092	PVC Thorax-drain Ch.32, straight, with adapter	10
401072	PVC Thorax-drain Ch. 32, straight, trocar compliant	10
401062	PVC Thorax-drain Ch.32, curved, with adapter	10
401057	PVC Thorax-drain Ch. 32, curved, trocar compliant	10

Ecolab Healthcare
Clean. Safe. Healthy.

Microtek Medical B.V.
A Division of Ecolab
Hekkehorst 24
NL-7207 BN Zutphen
Tel.: +31 (0) 575 599200

Ecolab Europe GmbH
Richtistr. 7
CH-8304 Wallisellen
www.ecolab.com



Drains

Closed Drainage Systems



Part.no.	Description	Qty/cs
Consisting of Silicon Drainage, x-ray positive with 5 eyes and 350 ml collection bag with tap, single packed, double peel-pack, sterile		
408510	Closed Drainage System Ch. 10	20
408512	Closed Drainage System Ch. 12	20
408515	Closed Drainage System Ch. 15	20
408518	Closed Drainage System Ch. 18	20
408521	Closed Drainage System Ch. 21	20
408524	Closed Drainage System Ch. 24	20
408527	Closed Drainage System Ch. 27	20
408530	Closed Drainage System Ch. 30	20
408536	Closed Drainage System Ch. 36	20

Drainage Systems Dual



Part.no.	Description	Qty/cs
Consisting of Silicon Drainage, (50 cm), x-ray positive with 5 eyes and 2000 ml collection bag with tap, double peel-pack, sterile		
408612	Drainage Systems Dual Ch. 12	20
408615	Drainage Systems Dual Ch. 15	20
408618	Drainage Systems Dual Ch. 18	20
408621	Drainage Systems Dual Ch. 21	20
408624	Drainage Systems Dual Ch. 24	20
408627	Drainage Systems Dual Ch. 27	20
408630	Drainage Systems Dual Ch. 30	20
408636	Drainage Systems Dual Ch. 36	20

Ecolab Healthcare
Clean. Safe. Healthy.

Microtek Medical B.V.
A Division of Ecolab
Hekkehorst 24
NL-7207 BN Zutphen
Tel.: +31 (0) 575 599200

Ecolab Europe GmbH
Richtistr. 7
CH-8304 Wallisellen
www.ecolab.com

