


Silicon Capillary-Drains

Thanks to the grooved internal wall, silicon capillary drains have an excellent capillary effect. These drains are generally used in general surgery and gynecology

Product description	Width	Article no.
Length 30 cm	6 mm	42 14 06
With x-ray	8 mm	42 14 08
	10 mm	42 14 10
	12 mm	42 14 12
Ditto, length 50 cm	12 mm	42 16 12



Silicon Thorax-Drains

Silicon thorax drains are specially designed for thorax and heart surgery.

The silicon thorax drains can be connected perfectly to Microtek Medical's thoracic drainage systems.

Product description				Article no.
With x-ray contrast				
Charr.	Length in cm	Perf. cm	Eyes	
18	50	12	8	40 11 18
24	50	17	10	40 11 24
27	60	20	10	40 11 27
30	65	24	12	40 11 30
36	70	27	14	40 11 36



Drainage systems

Closed



Contains Silicon-Drainage, x-ray compatible with 5 eyes and 350 ml collection bag with tap, double peel-pack, sterile


Product description	Article no.	
French/Charriere 10		40 85 10
	12	40 85 12
	15	40 85 15
	18	40 85 18
	21	40 85 21
	24	40 85 24
	27	40 85 27
	30	40 85 30
	36	40 85 36

Dual



Contains Silicon-Drainage, (50 cm), x-ray compatible with 5 eyes and 1.500 ml collection bag with tap, double peel-pack, sterile

Product description	Article no.	
French/Charriere 12		40 86 12
	15	40 86 15
	18	40 86 18
	21	40 86 21
	24	40 86 24
	27	40 86 27
	30	40 86 30
	36	40 86 36



Silicon Drains

Made from 100% medical silicon

Microtek Medical is your trusted source of innovative product solution for risk reduction and protection

Silicon Drains

Silicon Thoracic-Drains

Silicon Capillary-Drains

Silicon T-Drains

Silicon Flat-Drains




MICROTEK
MEDICAL B.V.
P.O. Box 234 7200 AE Zutphen The Netherlands
Hekkehorst 24 7207 BN Zutphen The Netherlands
Tel. +31 (0)575 599 200 Fax +31 (0)575 599 299

www.microtekmed.com

www.microtekmed.com



Silicon Drains made from 100% medical silicon



- smooth surface
- flexible
- moisture repellent
- excellent tissue and organ tolerance for longer periods of time
- excellent capillary effect due to the grooved internal wall
- rounded eyes and openings
- ETO sterilized

Microtek Medical's silicon drains are manufactured with the greatest care and are made of 100% medical silicon.

As a result, they are completely neutral, do not cause any tissue reactions and can be tolerated in organs and tissues for long periods of time.

Thanks to the moist-repellent properties of the silicon drains, tissue does not stick to the drain wall. Microtek Medical's silicon drains are individually wrapped in double-layer packaging and sterilized by gamma radiation.

These very flexible drains, which are not affected by temperature, are ready for use.

Silicon Drains

Silicon Drains for general Surgical use f.i. abdominal surgery


Product description	French/Charriere	Article no.
Length 50 cm, 5 lateral	8	40 51 08
holes with x-ray striping	10	40 51 10
	12	40 51 12
	15	40 51 15
	18	40 51 18
	21	40 51 21
	24	40 51 24
	27	40 51 27
	30	40 51 30
	36	40 51 36



Silicon T-Drains

For f.i. gall way drainage

Product description		Article no.
Silicon T-drains	3 mm	60 40 03
Silicon T-drains	4 mm	60 40 04
Silicon T-drains	5 mm	60 40 05
Silicon T-drains	6 mm	60 40 06
Silicon T-drains	7 mm	60 40 07
Silicon T-drains	8 mm	60 40 08



Flat drains

Flat drains type "Jackson Pratt" of medical silicon, x-ray positive, double-layer peel-pack, sterile

Product description			Article no.
Width	Length	Eyes	
4	65 cm	1-row	85 12 04
7	125 cm	2-rows	85 12 07
10	105 cm	2-rows	85 12 10
Suction balloon			
100cc	Silicon Bulb		85 11 20

