

**In USA:**

Microtek Medical, Inc.  
512 Lehmberg Road  
Columbus, MS 39702 USA  
Phone: (662) 327-1863  
Toll Free Phone: (800) 844-0988  
Inside Sales: (800) 866-8795  
Customer Service: (800) 824-3027  
Fax: (662) 327-5921

**In Europe and Middle East:**

Microtek Medical, B.V.  
Gerritsenweg 7  
7202 BP Zutphen  
The Netherlands  
Tel: +31 575 599 200  
Fax: +31 575 599 299  
UK Sales Office:  
Tel: +44 1782 561911  
Fax: +44 1782 562177

**All other countries:**

Microtek International Customer Service  
Phone: +1 (662) 244-3192  
U.S. Fax: +1 (662) 244-3198

[www.microtekmed.com](http://www.microtekmed.com)

© 2005 Microtek Medical, Inc. ISODrape® and Lingeman® are registered trademarks  
of Microtek Medical, Inc. Microban® is a registered trademark of Microban Products Company, Inc.



LIT013260



**Patient Drapes. Innovations  
In Protection With Features  
That Enhance Procedures.**

The Proven Leader In Prevention And Protection, Microtek Offers Advanced Patient Drapes For Risk Reduction And Performance Advantages For Urology and Obstetrics-Gynecology.



# Microtek. The trusted source of innovative patient draping for risk reduction and protection.



For more than two decades, Microtek Medical has been the trusted source of innovative patient draping solutions for risk reduction and protection. Microtek patient draping reduces the risks of infectious disease transmission; protects health-care professionals in every type of urological, obstetric and gynecological environment; and actually helps improve clinical outcomes by maximizing the protection available for both patients and medical staff during procedures.

There is no better example of this dedicated expertise than our complete range of patient drapes specially designed for urology and obstetrics-gynecology.

## Innovative thinking makes the difference in patient draping.

In urology, obstetrics and gynecology, innovative patient draping solutions from Microtek have significantly advanced the level of risk reduction and protection. This is not just

about producing drapes that meet needs. It's about designing more intelligent draping solutions that function effectively during every procedure. It's the difference made by extra features and thoughtful advances that can enhance the clinical environment and help improve outcomes.

## The special needs of urology and obstetrics-gynecology require something extra behind every patient drape: Microtek.

You already know Microtek as the long established leader in drapes for a diverse range of applications. The same thinking has been applied to patient draping designed exclusively for urology and obstetrics-gynecology. Where there are special



**"Innovation has long been a hallmark of the Microtek name. With every passing year, we can depend on new levels of ease, convenience and safety in every patient drape."**

needs from fluid control to higher levels of flexibility and functionality for select procedures, Microtek has created specialty specific patient drapes with added features to assure ease of use and exceed your every expectation.







## **Why draping solutions should be designed for patients and medical professionals.**

Each day, both healthcare professionals and patients come face to face with the serious risks of contamination and infection in urological and obstetric-gynecological environments. It's impossible to know every potential risk people may face in every clinical situation. It is, however, possible to protect every clinical environment. With a comprehensive line of advanced draping solutions specially designed for the procedures you perform on your patients, Microtek offers the key to peace of mind when it matters most.











**Order your Microtek patient drapes from the selection on the following pages, or visit [www.microtekmed.com](http://www.microtekmed.com) to see our complete online catalog.**

# C-Section Drapes

	Product #	Description	Qty/cs
	ISO6698PP	C-Section Drape with Viewing Window, 100" X 124" (254 X 315 cm) w/armboard covers, incise adhesive area 14" X 12" (36 X 32 cm) w/7.5" X 7.75" (19 X 20 cm) aperture, 4 cord holders, PerfectPouch™ fluid collection pouch, and viewing window 8" X 10" (20 X 25.4 cm)	5
	3460	Seattle C-Section Drape with Viewing Window, 100" X 112" (254 X 285 cm) w/armboard covers, incise adhesive area 14" X 12" (36 X 32 cm) w/7.5" X 7.75" (19 X 20 cm) aperture, 2 cord holders, U-shaped fluid collection pouch 45" X 30" (114 X 78 cm) w/2 drainage ports, and viewing window 8" X 10" (20 X 25.4 cm)	5
	ISO6697	C-Section Drape 100" X 124" (254 X 315 cm) w/armboard covers, square aperture 14" X 14" (36 X 36 cm) w/Microban antimicrobial incise adhesive, and fluid collection pouch 34" X 30" (86 X 76 cm) w/2 tube holders.	5
	3542	C-Birth Drape, 100" X 124" (254 X 315 cm) w/armboard covers, incise adhesive area 14" X 12" (36 X 32 cm) w/7.25" X 7.75" (18 X 20 cm) aperture, 2 cord holders, U-shaped fluid collection pouch 45" X 30" (114 X 78 cm) w/2 drainage ports. Non-stocking.	10
	3669	C-Birth Drape, 77" X 120" (196 X 305 cm) w/wings, square aperture 9" X 9" (23 X 23 cm) w/adhesive, fluid collection pouch 34" X 29" (86 X 74 cm) w/2 cord holders. Non-stocking.	10
	3736	C-Birth Drape, 77" X 120" (196 X 305 cm) w/wings, triangular aperture 10" X 11" (25 X 28 cm) w/incise adhesive, fluid collection pouch 34" X 29" (86 X 74 cm) w/3 cord holders. Non-stocking.	10

## Lingeman® Urology & GYN Specialty Products





	Product #	Description	Qty/cs
	1-0275	Lingeman Endo Bladder Suspension Drape 105" X 114" (267 X 290 cm) w/ 3" X 10" X 10" (8 X 25 X 25 cm) T-shaped aperture, leggings, armboard covers, 2-poly ties, 2 velcro cord holders, 14" X 10" (36 X 25 cm) flap w/ velcro, and fluid collection pouch 20" X 19" (51 X 48 cm) w/filter and drainage port. Also available non-sterile as product #1-0275N, packaged 25/case.	8

	1-0276	Lingeman Endo Bladder Suspension Pack contains: One Lingeman Endo Bladder Suspension Drape, One Suction Tubing 0.25" X 20' (0.6 X 610 cm) and One Back Table Cover/Outerwrap 50" X 90" (127 X 229 cm)	6
	1-0425PP	3 in 1 Drape 105" X 100" (267 X 254 cm) w/ 3" X 7" (8 X 18 cm) aperture w/ adhesive, leggings, 2 poly ties, PerfectPouch™ fluid collection pouch 24" X 18" (61 X 46 cm) w/ filter and drainage port.	10
	1-0425	3 in 1 Drape 105" X 100" (267 X 254cm) w/ 3" X 7" (8 X 18 cm) aperture w/ adhesive, leggings, 2 poly ties, fluid collection pouch 24" X 18" (61 X 46 cm) w/ filter and drainage port.	10
	1-0426	Lingeman 3 in 1 Pack Contains: One Lingeman 3 in 1 Drape, One suction tubing 0.25" X 12' (0.64 X 366 cm), and One Back Table Cover/Outerwrap 50" X 90" (127 X 229 cm)	6
	1-0657	Lingeman GYN Surgery Drape 107" X 130" (272 X 330 cm) w/10.5" X 10.5" (27 X 27 cm) abdominal aperture w/adhesive, 4" X 4" (10 X 10 cm) perineal aperture w/adhesive, clear anesthesia screen, armboard covers, leggings, 2-poly ties, cord holders, modesty flap, and fluid collection pouch 20" X 19" (51 X 48 cm) w/ filter and drainage port. 2 unattached dual laparoscopy instrument pockets.	6
	1-0658	Lingeman GYN Surgery Pack Contains: Two Absorbent Towels, 1 Mayo Stand Cover, 1 Under Buttocks Drape, 1 Suction Tubing 0.25" X 20' (0.64 X 610 cm), One Lingeman GYN Surgery Drape, 2 Dual Laparoscopy Instrument Pouches, and One Back Table Cover/Outerwrap 50" X 90" (127 X 229 cm)	4
	13-11131PP	Lingeman TUR Drape w/ Latex Free Finger Cot 105" X 110" (267 X 279 cm) w/ 3" X 7" (8 X 18 cm) aperture w/ adhesive, leggings, clear anesthesia screen, 2 poly ties, 1 latex free fingercot, and PerfectPouch™ fluid collection pouch 24" X 18" (61 X 46 cm) w/ filter and drainage port.	8




## How to order.

**For product information in the United States, call 800-886-8795 or to place an order, call 800-824-3027.**




For product information or to place an order in Europe and the Middle East, call our United Kingdom office at +44 1782 561911. In all other countries, call our U.S. office at +1 662-244-3192 or fax +1 662-244-3198. Or visit our website at [www.microtekmed.com](http://www.microtekmed.com) to review our online catalog, or contact your sales representative. Remember, the items featured in this catalog are just a small sample of the complete spectrum of Microtek Medical equipment draping solutions.

	13-11131	Lingeman TUR Drape w/ Latex Free Finger Cot 105" X 110" (267 X 279 cm) w/ 3" X 7" (8 X 18 cm) aperture w/ adhesive, leggings, clear anesthesia screen, 2 poly ties, 1 latex free fingercot, and fluid collection pouch 24" X 18" (61 X 46 cm) w/ filter and drainage port.	8
	13-11132	Lingeman TUR Pack Contains: Two Absorbent Towels, One Suction Tubing 0.25" X 20' (0.64 X 610 cm), One Lingeman TUR Drape, One Back Table Cover/Outerwrap 50" X 90" (127 X 229 cm)	6
	13-11134	Lingeman TUR Pack Contains: Two Absorbent Towels, One Suction Tubing 0.25" X 20' (0.64 X 610 cm), One 1' (30.5 cm) Tubing with Y Connectors, One Lingeman TUR Drape, and One Underbuttocks Pouch Drape*, One Back Table Cover/Outerwrap 50" X 90" (127 X 229 cm) *Covered by US Patent Number 6,314,958B1	6
	14-11173	Lingeman Lithotomy Endoscopy Drape with built-in Under Buttocks Drape, 95" X 101" (241 X 257 cm), w/clear anesthesia screen, leggings, 2 poly straps, w/3" X 7" (8 X 18 cm) aperture w/adhesive, built-in under buttocks pouch drape*, and fluid collection pouch w/filter and drainage port. *Covered by US Patent Number 6,314,958B1	8

## Microtek Patient Drapes For Urology and Obstetrics-Gynecology

	Product #	Description	Qty/cs
	ISO1071	ISODrape Urological Drape w/ latex-free fingercot 24" X 36" (61 X 91 cm) w/ 2.5" (6.35 cm) aperture.	10
	1-0815	Lingeman Percutaneous Procedure Drape 132" X 77" (335 X 196 cm) w/ 9" X 6" (23 X 15 cm) incise adhesive area with 1" (2.5 cm) aperture, 6 cord holders, fluid collection pouch 20" X 20" (51 X 51 cm) w/ filter and drainage port.	10
	1-0816	Lingeman Percutaneous Procedure Pack Contains: One C-Arm Drape 41" X 74" (104 X 188 cm), One Mayo Stand Cover 23" X 54" (58 X 137 cm), One Utility Table Cover 44" X 54" (112 X 137 cm), One Suction Tubing 0.25" X 12" (0.6 X 366 cm), One Lingeman Percutaneous Procedure Drape, One Back Table Cover/Outerwrap 50" X 90" (127 X 229 cm)	6

# Urology & OB/Gyn Pouches & Miscellaneous Drapes

	Product #	Description	Qty/cs
	ISO6658	ISODrape C-Birth Fluid Collection Pouch 29" X 34" (74 X 86 cm) Featuring 13" X 13" (33 X 33 cm) Microban Incise Adhesive area.	5
	2997	Peri/Ob Pouch Drape 39" X 46" (99 X 118 cm) w/underpad, full cuff and fluid collection pouch w/filter and drainage port.	20
	3885	Circumcision Drape 20" X 12" (51 x 30 cm) w/ 5/8" (1.58 cm) aperture, poly foam laminate.	40

## How to order.

**For product information in the United States, call 800-886-8795 or to place an order, call 800-824-3027.**

For product information or to place an order in Europe and the Middle East, call our United Kingdom office at +44 1782 561911. In all other countries, call our U.S. office at +1 662-244-3192 or fax +1 662-244-3198. Or visit our website at [www.microtekmed.com](http://www.microtekmed.com) to review our online catalog, or contact your sales representative. Remember, the items featured in this catalog are just a small sample of the complete spectrum of Microtek Medical equipment draping solutions.