



# Netis Disposable Cleaning Brushes



Our wide range of disposable cleaning brushes meets all your requirements for the perfect cleaning of your gastroscopes, colonoscopes, enteroscopes, nasofibrosopes and bronchoscopes.

### Features and Benefits:

- Easy to use: individually packed
- Easy to identify: color coded
- Easy to track: 3 traceability labels
- Reduced costs
- Optimum infection risk control
- Softness & efficiency



Description	Code-Number	Diameter Wire	Diameter Brush	Length	Color-Code
NETIS™ (with single head) for Colo/Gastroscopy	BNO-18230	1,7 mm	5,5 mm	230 cm	Green
NETIS™ dual (with double head) for Colo/Gastroscopy	DEC-18230	1,7 mm	5,5 mm	230 cm	Green
NETIS™ L for Enteroscopy	BNO-18280	1,7 mm	5,5 mm	280 cm	Green
NETIS™ C for Urology	BNO-18075	1,7 mm	5,5 mm	75 cm	Green
NETIS™ N for Naso-gastroscopy	BNO-17200	1,7 mm	3 mm	200 cm	White
NETIS™ P for Pneumology	BNO-17100	1,7 mm	3 mm	100 cm	White
NETIS™ W for Waterjet Kanal	BNO-10230	1,0 mm	2 mm	230 cm	Yellow
NETIS™ valv for Valve	BLP-10040		10 mm	11 cm	

### NETIS kit, Cleaning Kits

Description	Code-Number
BLP-10040 + BNO - 18230	KT-002
BLP-10040 + BNO-17200	KT-003
BLP-10040 + BNO-17100	KT-006
BLP-10040 + DEC-18230	KTE-002
BLP-1004- + DEC-17200	KTE-003
BLP-10040 + BNO-18230 + BNO-10230	KTP-005



P.O. Box 234, 7200 AE Zutphen  
Hekkehorst 24, 7207 BN Zutphen  
Tel. +31(0)575 599 200 Fax +31(0)575 599 299

[www.microtekmed.com](http://www.microtekmed.com)